


Developing arches with Transforce® appliances

Sagittal Expander



The Transforce® Arch Developers were created in consultation with Dr. William Clark, the inventor of the Clark Twin Block.

Transverse Expander



The Transforce offers aesthetics that adults prefer.

These appliances are virtually invisible in arch development, which appeal to adult patients. This system is also excellent for mixed dentition arch development.

The Transforce offers versatility.

All classes of malocclusions, as well as both upper and lower arches can be treated with Transforce appliances.

The Transforce is available in two styles.

The Transverse and the Sagittal.

The *Transverse* expands at the cuspids in both the upper and lower arches and is available in 2 upper and 2 lower sizes.

The *Sagittal* creates space by advancing the 3x3 mesially and is also available in 3 upper and 3 lower sizes.

The Twin force® style modules provide gentle biocompatible force.

Each expander contains nickel titanium springs generating approximately 200 grams of force. Both are easy to fit.

“Adjustments will be minimal to none with the spring-loaded Transforce system.”



Appliances that “snap-in” the first time.

800/457-0504 • www.johnsdental.com