



General Lab RY

PLEASE SEND
General Lab Rx
Ortho Rx
Ceph Rx
Shipping Boxes
Shipping Labels


423 South 13th St. • Terre Haute, IN 47807 • 800-457-0504
Local: (812) 232-6026 • Fax: (812) 234-4464 • Internet: www.johnsdental.com

CROWN & BRIDGE / PORCELAIN

NEXT APPOINTMENT

- Check Dies
Tryin Casting
Tryin Bisque Bake
Finish

Teeth # \_\_\_\_\_ High Noble
Teeth # \_\_\_\_\_ Noble
Teeth # \_\_\_\_\_ Base
Pontic # \_\_\_\_\_



- Porcelain Occlusal
Metal Occlusal
Metal Lingual Band
Circle metal design from metal design chart below.

MARYLAND BRIDGE

Teeth # \_\_\_\_\_ Abutments \_\_\_\_\_ Pontic \_\_\_\_\_

IMPLANTS

- Custom Abutments
Titanium
Zirconia
Attachments / Over Dentures
Locator
Hader Bar

ALL PORCELAIN OR COMPOSITE

- DJL Full Zirconia
Zirconium Oxide Substructures - Crown & Bridges
CeraC InLab Ditero
Porcelain Laminate Veneer
Demax


FULL CAST

Teeth # \_\_\_\_\_ High Noble
Teeth # \_\_\_\_\_ Noble
Teeth # \_\_\_\_\_ Base
Teeth # \_\_\_\_\_ Yellow
Teeth # \_\_\_\_\_ White


SURVEY FOR PARTIALS

- Survey for Partial
Saddle-Lock
Contour-to-Existing-Clasp (Send Partial)
2nd Nature Temporaries
Monodont Teeth

SPECIAL TECHNIQUES



SPECIAL INSTRUCTIONS



Porcelain Margin

www.johnsdental.com

LAB USE:
Please do not count Fridays, Saturdays, Sundays, holidays or days in transit as production days.

DATE SENT: / /

DATE WANTED: / /

DR. NAME: Last Name First Name

STREET ADDRESS:

CITY: STATE: ZIP:

DR. EMAIL: Last Name First Name

PATIENT NAME: Last Name First Name

AGE: I am a new customer

Male Female

Please contact me on this case

PARTIAL DENTURE

NEXT APPOINTMENT

- Upper Lower
Metal All-Acrylic (check one)
Framework Only: Tryin
Framework: Bite Rim
Framework: Set up with teeth
Flexible Partial: Set up with teeth
Flexible Partial: Set up and finish
Duraflex Clearnet
Valplast Acetal Resin
Alcast Veneers
Bite Restorer: Metal Occl.
Bite Restorer: Composite Occl.
Bite Restorer: Clear Acrylic Occl.
Bite Rim Finish (check one)
Posterior coverage only Antipost coverage

CLASPS

- As survey indicates
Akers
Bent Wire
RPI
Roach
Saddle-Lock
Tooth-Colored
Clearnet

CONNECTORS

- UPPER
Horseshoe
Palatal Bar
Double Palatal Bar
LOWER
Lingual Bar
Lingual Plate
Double Lingual Bar

MATERIAL
Vitalium 2000 (FDS Standard)
Gold
Design and Estimate Only (no charge)

TYPE OF CASE

- FULL PARTIAL IMMEDIATE
Upper Lower
Upper Lower

FULL DENTURE

Tooth Shade (Use tooth diagram if necessary)
Tooth Mould
Vitrapan
Set up: Vigorous Delicate

ESTHETICS

Posterior Tooth Type
Economy Teeth

SPECIAL TECHNIQUES

- Next Appointment ALL ACRYLIC
Custom Tray
Bite Rims
Control Rim Anterior Setup
Tracer
Setup Tryin
Surgical Stint

REPAIR, RELINE OR REPRODUCE

- Upper Lower
Soft Reline
Hard Reline
Repair
Reproduce

TMI & SLEEP APNEA

- Material
Ludacore 199 (Gum-like style Jit. weight)
ETHIC (Moderate)
Other

- Upper Lower
Thermo-Plastic (Clear Shell)
Hard Nightguard
Baker Deprogrammer
Sagittal
Pivotal: Standard
Special: Boiler Distal Push
Transitional: Myotonic Neuromuscular
Pull Forward: Myotonic Neuromuscular
Stuck Dentur (Hinge)
Farrar (w/ Anterior Plug)
Bryon Kemp
Clear Occlusal

- PRO-FORM SOFT MOUTHGUARDS
Upper Lower
Soft Nightguard
Standard Mouthguard
PowerBite Mouthguard
Heimert Strap (Optional)

- BLEACHING TRAYS
Upper Lower
Hard/Clear
Soft/Clear

- FINISH
Stipple
Smooth
Rugae
Palatal Relief
Posidam
Denture I.D.

- SLEEP APNEA
EVA
EVA (1st Step)
NAPA
SILENCER
SNOAR
TAP
LUCCO
HERBST
mRNA
OASYS

Accounts are due and payable upon receipt of monthly statement. All amounts not paid by the 23rd day of the month following the statement date are subject to a service charge on the unpaid balance at the rate of 2% per month. (2.4% per annum)
Accounts not paid within these credit terms will be subject to C.O.D. status.
Client bills, in full, the stated price of the goods, plus any service charges, plus all costs of collection including attorneys' fees, court costs & other reasonable expenses.

DRS. SIGNATURE:
REGISTRATION

LICENSE NO.

WHITE: Lab copy
YELLOW: Doctor's copy

Call us regarding biocompatible options.

www.johnsdental.com



800/457-0504





Cast Clasp Connected



Lingual Plate



Double Lingual Bar



Lingual Bar



Double Palatal Bar



Palatal Bar



Horseshoe



Saddle-Lock



RPI



Roach



Bent Wire



Akers

---- PARTIAL ----

**ClearMet** - Nylon-based, clear, flexible partial. The frame on the ClearMet Partial serves the same purpose as a metal frame partial, but with new attractive features. The ClearMet is made of a material that is monomer free, FDA approved, and with no known allergic reactions. The partial is nonporous, resulting in a product that has no odors and is bacteria and stain resistant.

**Framework** - Unless otherwise requested, all frameworks are constructed with Vitalium 2000 giving doctors and patients a thinner framework design with added strength required for predictable adjustments.

**Alcast Veneers** - These veneers are custom shaded and shaped composite veneers placed on an Alcast framework. (Please indicate tooth# and shade.)

**Bite Restorer: Metal Occlusal** - A metal framework is used to cover posterior teeth and open the bite

**Bite Restorer: Composite occlusal** - (standard) A metal framework is used to open the bite, covering posterior as well as optional anterior teeth. The occlusal surface consists of a composite material. (choose shade) It is strongly recommended that a frame with bite rim be tried in before finishing.

**Bite Restorer: Clear Acrylic Occlusal** - Clear acrylic can also be used to cover posterior teeth as a lower cost alternative until a more stable bite is established. A more esthetic composite can replace the clear acrylic when the bite is stabilized.

**DuraFlex Partial** - DuraFlex flexible partials offer superb esthetics, hygiene and fit.

**Tooth-Colored Acetal Resin** - is an engineering thermoplastic used in precision parts requiring high stiffness, low friction and excellent stability. We have used them in partial frames and tooth-colored clasps when requested.

**Attachments**

Attachment systems recommended by Johns Dental Laboratories include Locator, ERA and Hader Bar.

**Implants**

Johns Dental is certified with all the major implant companies.

---- DENTURE ----

**AED (Avoid Embarrassment Denture)** - an economical spare denture made with a new denture.

**Twin Denture** - A simple system for building a transitional, immediate denture. The original setup is duplicated for immediate use (AED) and wax try-in is retained for processing after healing is complete.

**Virginia Partial** - This over-denture, with soft elastometric band(s), is made of Luci-Soft material that gently clasps and cushions select abutment teeth. The Virginia Partial is also a great way to cover visible roots common with receding tissue.

**Repairs, Relines, or Reproduce**

Most hard relines and repairs, received by 9:30 a.m., can be shipped in the same day. Soft relines must also be received in the morning but require overnight processing. A reproduce replaces an existing denture base, but retains the teeth and set up. It requires two days in the laboratory.

**Luci-Soft** - This heat-cured silicone based soft-liner offers long lasting softness, stronger bonding and stain-resistance.

---- TMJ SPLINTS ----

**Bryan Ramp** - The Bryan Ramp Therapeutic Night Guard is a major step forward in the treatment of chronic muscle tension from nocturnal clenching and grinding. Its FDA approved non-allergenic materials allow it to be thin yet strong.

**Clear Splint** - This methyl methacrylate free, non-yellowing material is an excellent material for splints that need to be flexible enough for a great fit, yet tough enough to withstand the rigors of bruxism. The thermoplastic material used in our Clear Splints softens in warm water.

**Denar/Witzig** - This mandibular splint is designed with full occlusal coverage and deep occlusal indentation that guides the patient into the proper bite established by the Denar/Witzig Articulator. A Night time Corrector may be used during evening wear.

**Farrar** - maxillary pull forward splint with incline plane in anterior section - posterior contact is between the buccal cusp of the lower teeth and the upper acrylic occlusal pad

**Gelb** - mandibular splint with acrylic coverage over posteriors with upper lingual cusps contacting a lower flat or indented occlusal pad (Metal Lingual Bar)

**Hard Nightguard** - Hard acrylic nightguard made on either upper or lower to provide relief for the bruxing and grinding. Send upper or lower models if indexing is desired.

**Myotronic** - mandibular splint with clear acrylic coverage over all teeth, built to a Myoprint bite and is made to look like natural dentition - mounted on a Terminus articulator

**Stack (Pull Forward)** - mandibular splint used to hold mandible forward and open for cases that click or were closed-locked and have been treated to a click

**Olmos Series** - Johns Dental has become a certified lab for the T&S Therapy Centre and offers a full line of Olmos Appliances. Please call to find out how you can become a certified provider.

**Pro-Form Soft Mouthguards**

**Soft: Mouthguard** - provides special protection for athletes in contact sports and comes in clear unless color is specified. Most basic colors and combinations are available. See website for color choices: [www.johnsdental.com](http://www.johnsdental.com)

**Soft: Nightguard** - soft material to provide relief from bruxing and grinding

**PowerBite Mouthguard** - This premium mouthguard combines protection with increased athletic strength & endurance.

---- SLEEP APNEA ----

**EMA (Elastic Mandibular Advancement)** - This sleep apnea appliance is formed from a dental impression and resembles and fits like a retainer. The use of elastic bands increases movement of the mandible, reduces the risk of TMJ discomfort and allows for individualized treatment. The major design innovation includes the regulation of the amount of advancement and freedom of mandibular movement through different bands.

**HERBST Sleep Appliance** - offers a telescopic mechanism which advances the mandible into a forward position allowing for an improved airway.

**LUCO** - the LUCO Hybrid OSA is unique in that it was originally designed and used for the treatment of posttraumatic TMD cases as a night stabilizing appliance to facilitate healing. It was not until a sleep specialist discovered that it was also very effective in treating sleep apnea that it was developed for both uses.

**mRNA** - Johns Dental is a certified lab for mRNA and DNA appliances. Please call for more information.

**NAIPA (Nocturnal Airway Patency Appliance)** - The NAPA reduces OSA and snoring by increasing tonicity of tongue, inhibiting movement of muscles of mastication, stabilizing mandibular anchorage of geniohyoid muscle and displacing the tongue anteriorly.

**OSA3YS** - This system functions by repositioning the mandible in addition to acting as a nasal dilator for reduction of nasal resistance and improved nasal breathing.

**Silencer** - The Silencer is a 2-component, articulated snoring & sleep apnea appliance that advances the mandible through a controlled range of protrusive movements. The Silencer uses a commercially pure titanium hinge to provide incremental advancement of the mandible.

**Snoozer** - The Snoozer appliance is a one-piece orthodontic appliance that moves the lower jaw forward along with the tongue, which opens the airway and facilitates normal breathing. Johns Dental is the sole and exclusive laboratory in the U.S. licensed to fabricate this appliance.

**TAP** - The TAP sleep apnea appliance maintains a clear airway to reduce snoring and improve breathing by holding the lower jaw in a forward position. (Write on Rx if the TAP case is for a Medicare patient. Appliance designs may be affected.)

---- CROWN AND BRIDGE ----

**Laminate Veneer** - full porcelain veneer covering the labial surface of anterior teeth with no metal substructure - most often prescribed for cosmetic enhancement

**Maryland Bridge** - semi-precious bridgework retained by lingual surfaces on abutment teeth - lingual surfaces are silicoated for the best possible retention

**Monodont** - Monodont has a variety of applications that involve mid-term fixed prostheses for both anterior and posterior applications. Monodont is structurally very durable, utilizing the simple engineering principle of "re-bar."

**Zirconium Oxide Substructures** - offers a stronger coping due to the densely-sintered zirconium oxide designed to withstand posterior demands, which provides doctors with another choice for anterior applications.

**JDL Full Zirconia** - A CAD/CAM procedure for the fabrication of all-ceramic crowns and bridges for anterior and posterior applications.

**Second Nature Temporaries** - Mid-term temporary crown/bridge - heat and pressure-cured method (durable & strong) - light-cured glaze resists stains

**Biocompatible** - Johns Dental offers several materials and techniques for those patients who require specific allergy solutions. Call us about specific case questions.