

Johns Dental Laboratories

Johns Dental Lab is a full service, family owned and operated dental laboratory that began over 75 years ago. We take pride in our finely crafted custom-made products, and our promise to doctors is that our appliances and prosthetics will fit your model.



Why is the ClearMet partial superior to other partials?

- ◆ Extremely Durable
- ◆ Easily Adjustable
- ◆ Chairside Relines
- ◆ Stain Resistant-no odor
- ◆ Invisible
- ◆ Tooth Bearing Clasp
- ◆ Add Additional Teeth-quickly and easily
- ◆ Amazing Comfort
- ◆ FDA approved

When we say we have a partial that disappears in the mouth, we aren't kidding!

Call us today to schedule a pick up!

Johns Dental Laboratories

423 South 13th Street
Terre Haute, IN 47807

Phone: 800-457-0504
Fax: 812-234-4464
www.johnsdental.com

ClearMet

The only clear choice for partials.



Appliances and Prosthetics That Fit Right the First Time.



ClearMet

The ClearMet Partial is the most functional and superior partial denture in the market today. The frame on the ClearMet Partial serves the same purpose as the metal framed partial, but with new attractive features that do not exist on any other partial.

ClearMet Partials are made of a material that is monomer free, FDA approved, and with no known allergic reactions. The partial is nonporous, resulting in a product that has no odors and is bacteria and stain resistant.

VS. Other Flexible Partials

Unlike conventional, flexible partials, the ClearMet Partial is produced using a two step procedure.

First, the frame is made through injection of the materials. The saddles are then processed onto the frame, which was designed on a surveyor to produce the most accurate final product as specified by the Provider.

We have the ability to create occlusal bite openers using the ClearMet material that cannot be achieved with traditional pink material. ClearMet Partials are easily repaired in the unlikely event of a break and have the ability to add teeth quickly and easily.

"I couldn't afford implants and I had no other choice but to get a partial denture for my missing teeth. I have had a metal partial for a few years now so I immediately expressed to my dentist my fear of smiling and seeing a metal object in my mouth. He put my nerves at ease explaining to me about the ClearMet partial and how they are invisible in the patient's mouth. I have had my new partial for almost a year now and I can't believe how nice it is to not have to think before smiling. I would recommend the ClearMet partial to anyone needing a partial denture." Shirley W.

Our Goals

Our goal for the patient is not only to improve smiles, but more importantly, improve lives. Not only will your patients be more confident in their appearance, but they will eat better, impacting their overall health.

Our goal for you, the Doctor, is to provide you with a superior partial denture, crafted to your patient's personal treatment needs.



Unlike other partials, ClearMet Partials are great for temporary implant restoration due to tooth bearing clasps. The tooth bearing clasps allow the tissue to be covered without coming in contact with the implant.