


EMA Strap lengths

SOFT **MEDIUM** **FIRM** **EXTRA FIRM**
 White Yellow Blue Clear

NEW STRAP CONVERSION CHART	
OLD CHART	NEW CHART
#1	21mm
#2	19mm
#3	17mm
#4	15mm*
#5	14mm

*The #4 straps prior to 2011 were 15.5 mm



SAMPLE PACK

A sample pack consisting of two pairs of white, yellow and blue of sizes 14, 17 and 21 will be included with every case.

The EMA appliance allows several simple variations to customize your treatment plan. These color-coded elastic straps provide varying degrees of mandibular advancement. Your patient will also be able to maintain unsurpassed lateral movement and overall TMJ comfort.

Call Johns Dental at 800/457-0504 to order additional straps.

www.johnsdental.com

(Johns Dental places buttons 25 mm apart.)

